

Human Apolipoprotein E (APOE) ELISA Kit

Catalog No: #EK20865

Package Size: #EK20865-1 48T #EK20865-2 96T

Description

Product Name	Human Apolipoprotein E (APOE) ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Apo-E; AD2; Apoprotein; Alzheimer Disease 2(E4-Associated,Late Onset

Application Details

Detect Range:23.4-1500ng/mL

Sensitivity:8.47ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.