

# Human Deoxyribonuclease I Like Protein 3 (DNASE1L3) ELISA Kit

Catalog No: #EK20866

Package Size: #EK20866-1 48T #EK20866-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                               |
|--------------------|---------------------------------------------------------------|
| Product Name       | Human Deoxyribonuclease I Like Protein 3 (DNASE1L3) ELISA Kit |
| Applications       | ELISA                                                         |
| Species Reactivity | Human                                                         |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.065ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.