

Mouse Aldehyde dehydrogenase family 1 member A3(ALDH1A3) ELISA Kit

Catalog No: #EK20873

Package Size: #EK20873-1 48T #EK20873-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aldehyde dehydrogenase family 1 member A3(ALDH1A3) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.32ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.