

Human Adhesion G protein-coupled receptor G3(ADGRG3) ELISA Kit

Catalog No: #EK20875

Package Size: #EK20875-1 48T #EK20875-2 96T

Description

Product Name	Human Adhesion G protein-coupled receptor G3(ADGRG3) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:15.6-1000pg/m

Sensitivity:5.8pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.