

Mouse Growth arrest and DNA damage-inducible protein GADD34(GADD34) ELISA Kit

Catalog No: #EK20879

Package Size: #EK20879-1 48T #EK20879-2 96T

Description

Product Name	Mouse Growth arrest and DNA damage-inducible protein GADD34(GADD34) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:11.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.