

Rat Pyruvate Dehydrogenase Beta (PDHb) ELISA Kit

Catalog No: #EK20880

Package Size: #EK20880-1 48T #EK20880-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Pyruvate Dehydrogenase Beta (PDHb) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.625-40ng/mL

Sensitivity:0.292ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.