

Human Transmembrane emp24 domain-containing protein 10(TMED10) ELISA Kit

Catalog No: #EK20881

Package Size: #EK20881-1 48T #EK20881-2 96T



Description

Product Name	Human Transmembrane emp24 domain-containing protein 10(TMED10) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.124ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.