

Human Arachidonate 12-lipoxygenase, 12R-type(ALOX12B) ELISA Kit

Catalog No: #EK20882

Package Size: #EK20882-1 48T #EK20882-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Arachidonate 12-lipoxygenase, 12R-type(ALOX12B) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range: 0.156-10 ng/mL

Sensitivity: 0.067 ng/mL

Sample Type: Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 µL

Assay Time: 1-4.5 h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.