

Rat Cyclic AMP-responsive element-binding protein 1(CREB1) ELISA Kit

Catalog No: #EK20889

Package Size: #EK20889-1 48T #EK20889-2 96T

Description

Product Name	Rat Cyclic AMP-responsive element-binding protein 1(CREB1) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:11.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.