

Mouse Transcription factor SOX-9 (SOX9) ELISA Kit

Catalog No: #EK20892

Package Size: #EK20892-1 48T #EK20892-2 96T



Description

Product Name	Mouse Transcription factor SOX-9 (SOX9) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.103ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.