

Rat Somatostatin (SST) ELISA Kit

Catalog No: #EK20897

Package Size: #EK20897-1 48T #EK20897-2 96T

Description

Product Name	Rat Somatostatin (SST) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:7.81-500pg/mL

Sensitivity:2.16pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.