

# Human Acetylcholinesterase collagenic tail peptide (COLQ) ELISA Kit

Catalog No: #EK20953

Package Size: #EK20953-1 96T #EK20953-2 48T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                     |
|--------------------|---------------------------------------------------------------------|
| Product Name       | Human Acetylcholinesterase collagenic tail peptide (COLQ) ELISA Kit |
| Applications       | ELISA                                                               |
| Species Reactivity | Human                                                               |

## Application Details

Detect Range:0.156ng/mL-10ng/mL Sensitivity:0.078 ng/mL Sample Type:Serum, Plasma, Other biological fluids Sample Volume: 1-200  $\mu$ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.