

Mouse Free fatty acids (FFA) ELISA Kit

Catalog No: #EK20954

Package Size: #EK20954-1 48T #EK20954-2 96T



Description

Product Name	Mouse Free fatty acids (FFA) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78.1 pg/mL - 5000 pg/mL Sensitivity:39 pg/mL Sample Type:Serum, Plasma, Other biological fluids Sample Volume: 1-200 μ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.