

Mouse Interferon alpha-5 (IFNA5) ELISA Kit

Catalog No: #EK20955

Package Size: #EK20955-1 48T #EK20955-2 96T



Description

Product Name	Mouse Interferon alpha-5 (IFNA5) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	RP11-380P16.5; INFA5;

Application Details

Detect Range:31.25 pg/mL - 2000 pg/mL Sensitivity:15.6 pg/mL Sample Type:Serum, Plasma, Other biological fluids Sample Volume: 1-200 μ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.