

Mouse Peptidoglycan Recognition Protein 1
(PGLYRP1) ELISA Kit

Catalog No: #EK20962

Package Size: #EK20962-1 48T #EK20962-2 96T



Description

Product Name	Mouse Peptidoglycan Recognition Protein 1 (PGLYRP1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL Sensitivity:5.9pg/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.