

Human Isthmin-1 (ISM1) ELISA Kit

Catalog No: #EK20963

Package Size: #EK20963-1 48T #EK20963-2 96T



Description

Product Name	Human Isthmin-1 (ISM1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL Sensitivity:28pg/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.