

Cat Interleukin 2 (IL-2) ELISA Kit

Catalog No: #EK20966

Package Size: #EK20966-1 48T #EK20966-2 96T



Description

Product Name	Cat Interleukin 2 (IL-2) ELISA Kit
Applications	ELISA
Species Reactivity	Cat
Other Names	IL-2; TCGF; lymphokine; T cell growth factor[aldesleukin]interleukin-2involved in regulation of T-cell clonal expansion

Application Details

Detect Range:15.6-1000 pg/mL Sensitivity:7.8 pg/mL Sample Type:Serum, Plasma, Other biological fluids Sample Volume:1-200 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.