

Human Secretory Immunoglobulin A (sIgA)ELISA Kit

Catalog No: #EK20967

Package Size: #EK20967-1 48T #EK20967-2 96T



Description

Product Name	Human Secretory Immunoglobulin A (sIgA)ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL Sensitivity:0.146ng Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.