

Human FRAS1-related extracellular matrix protein
2oO FREM2oO ELISA Kit

Catalog No: #EK20969

Package Size: #EK20969-1 48T #EK20969-2 96T



Description

Product Name	Human FRAS1-related extracellular matrix protein 2oO FREM2oO ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.053ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.