

Mouse Inhibitor of nuclear factor kappa-B kinase subunit alpha (IKKA) ELISA Kit

Catalog No: #EK20970

Package Size: #EK20970-1 48T #EK20970-2 96T



Description

|                    |                                                                                 |
|--------------------|---------------------------------------------------------------------------------|
| Product Name       | Mouse Inhibitor of nuclear factor kappa-B kinase subunit alpha (IKKA) ELISA Kit |
| Applications       | ELISA                                                                           |
| Species Reactivity | Mouse                                                                           |

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.057ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50  $\mu$ L Assay Time:1-4.5h  
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.