

Mouse Inhibitor of nuclear factor kappa-B kinase subunit beta (IKKB) ELISA Kit

Catalog No: #EK20971

Package Size: #EK20971-1 48T #EK20971-2 96T



Description

Product Name	Mouse Inhibitor of nuclear factor kappa-B kinase subunit beta (IKKB) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.054ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.