

Human Gamma-aminobutyric acid
receptor-associated protein(GBRAP) ELISA Kit

Catalog No: #EK20973

Package Size: #EK20973-1 48T #EK20973-2 96T



Description

Product Name	Human Gamma-aminobutyric acid receptor-associated protein(GBRAP) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL Sensitivity:0.118ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.