

Human von Willebrand factor C and EGF domain-containing protein(VWCE)ELISA Kit

Catalog No: #EK20975

Package Size: #EK20975-1 48T #EK20975-2 96T

Description

Product Name	Human von Willebrand factor C and EGF domain-containing protein(VWCE)ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.137ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ OL

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.