

# Human Olfactomedin-2(OLFM2) ELISA Kit

Catalog No: #EK20978

Package Size: #EK20978-1 48T #EK20978-2 96T

## Description

|                    |                                       |
|--------------------|---------------------------------------|
| Product Name       | Human Olfactomedin-2(OLFM2) ELISA Kit |
| Applications       | ELISA                                 |
| Species Reactivity | Human                                 |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.067ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.