

Human Voltage-dependent anion-selective channel protein 3 (VDAC3) ELISA Kit

Catalog No: #EK20981

Package Size: #EK20981-1 48T #EK20981-2 96T

Description

Product Name	Human Voltage-dependent anion-selective channel protein 3 (VDAC3) ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	HD-VDAC3

Application Details

Detect Range:0.156 ng/mL - 10 ng/mL Sensitivity:0.1 ng/mL Sample Type:Serum, Plasma, Other biological fluids Sample Volume:1-200 μ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.