

Human Trophoblast glycoprotein(TPBG) ELISA Kit

Catalog No: #EK20982

Package Size: #EK20982-1 48T #EK20982-2 96T

Description

Product Name	Human Trophoblast glycoprotein(TPBG) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.58pg/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.