

Human Ectonucleotide
pyrophosphatase/phosphodiesterase family member
3(ENPP3) ELISA Kit

Catalog No: #EK20984



Package Size: #EK20984-1 48T #EK20984-2 96T

Description

Product Name	Human Ectonucleotide pyrophosphatase/phosphodiesterase family member 3(ENPP3) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.625-40ng/mL Sensitivity:0.253ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.