

# Human Acylsphingosine deacylase 2(ASAH2) ELISA Kit

Catalog No: #EK20985

Package Size: #EK20985-1 48T #EK20985-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                    |
|--------------------|----------------------------------------------------|
| Product Name       | Human Acylsphingosine deacylase 2(ASAH2) ELISA Kit |
| Applications       | ELISA                                              |
| Species Reactivity | Human                                              |

## Application Details

Detect Range:0.625-40ng/mL

Sensitivity:0.225ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.