

Mouse R-Spondin 3 (RSPO3) ELISA Kit

Catalog No: #EK20986

Package Size: #EK20986-1 598 #EK20986-2 768



Description

Product Name	Mouse R-Spondin 3 (RSPO3) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.061ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.