

Human Four and a half LIM domains protein 2(FHL2) ELISA Kit

Catalog No: #EK20989

Package Size: #EK20989-1 48T #EK20989-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Four and a half LIM domains protein 2(FHL2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.064ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.