

Human Sorting nexin-5(SNX5) ELISA Kit

Catalog No: #EK20991

Package Size: #EK20991-1 48T #EK20991-2 96T

Description

Product Name	Human Sorting nexin-5(SNX5) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.