

Human Aldehyde Dehydrogenase Family 1 Member A3 (ALDH1A3) ELISA Kit

Catalog No: #EK20993

Package Size: #EK20993-1 48T #EK20993-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Aldehyde Dehydrogenase Family 1 Member A3 (ALDH1A3) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:31.2-2000pg/mL Sensitivity:11.6pg/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.