

Human Immunoglobulin heavy constant mu(IGHM)ELISA Kit

Catalog No: #EK20997

Package Size: #EK20997-1 48T #EK20997-2 96T

Description

Product Name	Human Immunoglobulin heavy constant mu(IGHM)ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.37ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ OL

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.