

Human Biogenesis of lysosome-related organelles complex 1 subunit 6(BLOC1S6) ELISA Kit

Catalog No: #EK20998

Package Size: #EK20998-1 48T #EK20998-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Biogenesis of lysosome-related organelles complex 1 subunit 6(BLOC1S6) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL Sensitivity:0.115ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.