

Human Endoribonuclease LACTB2(LACTB2) ELISA Kit

Catalog No: #EK20999

Package Size: #EK20999-1 48T #EK20999-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Endoribonuclease LACTB2(LACTB2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ OL

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.