

Human Cytochrome P450 26B1(CYP26B1) ELISA Kit

Catalog No: #EK21000

Package Size: #EK21000-1 48T #EK21000-2 96T



Description

Product Name	Human Cytochrome P450 26B1(CYP26B1) ELISA Kit
Applications	ELISA
Species Reactivity	Human
Crossing Reactivity	Detection Method:sandwich

Application Details

Detect Range:0.156-10ng/mL/br> Sensitivity:0.058ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.