

Human GNP NAT1 / Glucosamine 6-phosphate
N-acetyltransferase ELISA Kit

Catalog No: #EK21001

Package Size: #EK21001 96T



Description

Product Name	Human GNP NAT1 / Glucosamine 6-phosphate N-acetyltransferase ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Phosphoglucosamine acetylase,Phosphoglucosamine transacetylase,GNA1
Accession No.	Uniprot:Q96EK6
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.101ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.