

# Human DUOX1 / Dual oxidase 1 ELISA Kit

Catalog No: #EK21004

Package Size: #EK21004 96T

## Description

|                    |                                                                                   |
|--------------------|-----------------------------------------------------------------------------------|
| Product Name       | Human DUOX1 / Dual oxidase 1 ELISA Kit                                            |
| Applications       | ELISA                                                                             |
| Species Reactivity | Human                                                                             |
| Other Names        | Large NOX 1,Long NOX 1,NADPH thyroid oxidase 1,Thyroid oxidase 1,DUOX,LNOX1,THOX1 |
| Accession No.      | Uniprot:Q9NRD9                                                                    |
| Storage            | 6 months                                                                          |

## Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.066ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.