

Human VASH2 / TubulinyI-Tyr carboxypeptidase 2 ELISA Kit

Catalog No: #EK21006

Package Size: #EK21006 96T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human VASH2 / TubulinyI-Tyr carboxypeptidase 2 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	VASHL,Vasohibin-2,Vasohibin-like protein
Accession No.	Uniprot:Q86V25
Storage	6 months

Application Details

Detection Range:0.78-50ng/mLAssay Type:SandwichSensitivity:0.24ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.