

Human CCDC25 / Coiled-coil domain-containing protein 25 ELISA Kit

Catalog No: #EK21008

Package Size: #EK21008 96T



Description

Product Name	Human CCDC25 / Coiled-coil domain-containing protein 25 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:Q86WR0
Storage	6 months

Application Details

Detection Range:0.156-10ng/mLAssay Type:SandwichSensitivity:0.086ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.