

# Human ILF2 / Interleukin enhancer-binding factor 2 ELISA Kit

Catalog No: #EK21013

Package Size: #EK21013 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                              |
|--------------------|--------------------------------------------------------------|
| Product Name       | Human ILF2 / Interleukin enhancer-binding factor 2 ELISA Kit |
| Applications       | ELISA                                                        |
| Species Reactivity | Human                                                        |
| Other Names        | Nuclear factor of activated T-cells 45 kDa,NF45,PRO3063      |
| Accession No.      | Uniprot:Q12905                                               |
| Storage            | 6 months                                                     |

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.085ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.