

Human B9D1 / B9 domain-containing protein 1 ELISA Kit

Catalog No: #EK21014

Package Size: #EK21014 96T



Description

Product Name	Human B9D1 / B9 domain-containing protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	MKS1-related protein 1,MKSR1
Accession No.	Uniprot:Q9UPM9
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.085ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.