

Human ACMSD / 2-amino-3-carboxymuconate-6-semialdehyde decarboxylase ELISA Kit

Catalog No: #EK21017

Package Size: #EK21017 96T

Description

Product Name	Human ACMSD / 2-amino-3-carboxymuconate-6-semialdehyde decarboxylase ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Picolinate carboxylase
Accession No.	Uniprot:Q8TDX5
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.29ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.