

Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit

Catalog No: #EK21018

Package Size: #EK21018 96T

Description

Product Name	Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell recognition molecule Caspr4,CASPR4,KIAA1763
Accession No.	Uniprot:Q9C0A0
Storage	6 months

Application Details

Detection Range:0.156-10ng/mLAssay Type:SandwichSensitivity:0.056ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.