

Human CSRP3 / Cysteine and glycine-rich protein 3
ELISA Kit

Catalog No: #EK21020

Package Size: #EK21020 96T



Description

Product Name	Human CSRP3 / Cysteine and glycine-rich protein 3 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cardiac LIM protein,Cysteine-rich protein 3,CRP3,LIM domain protein, cardiac,Muscle LIM protein,CLP,MLP
Accession No.	Uniprot:P50461
Storage	6 months

Application Details

Detection Range:0.156-10ng/mLAssay Type:SandwichSensitivity:0.078ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.