

Human PEAK1 / Inactive tyrosine-protein kinase
PEAK1 ELISA Kit

Catalog No: #EK21024

Package Size: #EK21024 96T



Description

Product Name	Human PEAK1 / Inactive tyrosine-protein kinase PEAK1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	KIAA2002,Pseudopodium-enriched atypical kinase 1,Sugen kinase 269,Tyrosine-protein kinase SgK269
Accession No.	Uniprot:Q9H792
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:24pg/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.