

Human SLC6A1 / Sodium- and chloride-dependent GABA transporter 1 ELISA Kit

Catalog No: #EK21030

Package Size: #EK21030 96T

Description

Product Name	Human SLC6A1 / Sodium- and chloride-dependent GABA transporter 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	GAT-1, Solute carrier family 6 member 1, GABATR, GABT1, GAT1
Accession No.	Uniprot:P30531
Storage	6 months

Application Details

Detection Range: 78-5000 pg/mL Assay Type: Sandwich Sensitivity: 39 pg/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.