

Mouse Nlgn1 / Neuroligin-1 ELISA Kit

Catalog No: #EK21033

Package Size: #EK21033 96T



Description

Product Name	Mouse Nlgn1 / Neuroligin-1 ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Kiaa1070
Accession No.	Uniprot:Q99K10
Storage	6 months

Application Details

Detection Range:0.156-10ng/mLAssay Type:SandwichSensitivity:0.078ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.