

Human UBXN7 / UBX domain-containing protein 7 ELISA Kit

Catalog No: #EK21034

Package Size: #EK21034 96T

Description

Product Name	Human UBXN7 / UBX domain-containing protein 7 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	KIAA0794,UBXD7
Accession No.	Uniprot:O94888
Storage	6 months

Application Details

Detection Range:0.78-50ng/mLAssay Type:SandwichSensitivity:0.41ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.