

Human GKN1 / Gastrokeine-1 ELISA Kit

Catalog No: #EK21035

Package Size: #EK21035 96T



Description

Product Name	Human GKN1 / Gastrokeine-1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	18 kDa antrum mucosa protein,AMP-18,Protein CA11,AMP18,CA11,UNQ489/PRO1005
Accession No.	Uniprot:Q9NS71
Storage	6 months

Application Details

Detection Range:0.312-20ng/mLAssay Type:SandwichSensitivity:0.186ng/mL
--

Note: This product is for in vitro research use only and is not intended for use in humans or animals.